

Introducing the...

Coda® Canister Main Stage Filter

IVF Air Quality enters the 21st Century... for you and your embryos

The Coda® Canister air filters gives you more clean activated carbon and impregnated potassium permanganate than the competition.



Coda® Canister...

Developed for IVF Air Quality

- Maximum absorption
- More residence time
- More highly activated, clean carbon
- More impregnated potassium permanganate
- Greater removal of air contaminants

The Coda® Canister Main Stage Filter gives you 7 important reasons to use this filter over all others.

1. Superior absorption capacity compared to any other filter, per carbon pound.
2. Air tight seals for single and multiple pass systems.
3. Performs for up to 6 months removing the majority of VOCs and CACs.
4. Will maintain a higher air exchange rate throughout the life of the filter.
5. Removes the high amounts of Volatile Organic Compounds.
6. Removes the high amounts of Chemical Air Contaminants.
7. Economical & Efficient.



FDA 510(k) Cleared

<genX>
international




IVFonline, LLC
Leading Fertility Brands

The Coda® Canister Main Stage Filter

The Coda® Canister Main Stage Filter begins with a >95% absorption efficiency rate which is far above the 57% efficiency rate of conventional loose carbon filled V-cell filters.

Coda® Canisters, installed 5 in our custom frame, will adsorb almost 1 lb of VOCs and CACs, 2 times that of a V-cell filter while still maintaining the same efficiency rate.



The Coda® Canister will continually remove more harmful contaminants from the air in your IVF laboratory, without a reduction in air exchanges, protecting you and your specimens.

The Coda® Canister easily loads into the Coda Tower®, Coda® Lo-Boy and most HVAC systems.

Item	Frequency	Catalogue #
Coda® Canisters (5)	Six Months	CCST-005
Coda® Canister Frame	Reusable*	CCST-000
Coda® Canister Initial Kit Consists of: <ul style="list-style-type: none"> • 5 Coda® Canisters • 1 Coda® Canister Frame • 1 Coda® HEPA Filter • 2 Coda® PreFilters 	Initial	For Coda Tower®: CCTI-081
		For Coda® Lo-Boy: CCLI-082
Coda® Canister Semi-Annual Kit Consists of: <ul style="list-style-type: none"> • 5 Coda® Canisters • 2 Coda® PreFilters 	Semi-Annual	For Coda Tower®: CCTS-081
		For Coda® Lo-Boy: CCLS-082
Coda® Canister Annual Kit Consists of: <ul style="list-style-type: none"> • 5 Coda® Canisters • 1 Coda® HEPA Filter • 2 Coda® PreFilters 	Annual	For Coda Tower®: CCTA-081
		For Coda® Lo-Boy: CCLA-082

*The frame is reusable and you should order one at the initial purchase of the Coda® Canister System.

Worldwide Service at Your Fingertips

Call or Click



US/Canada: 800-720-6375
International: 519-826-5800

www.IVFOonline.com

US/Canada: 800-720-6375 • International: 519-826-5800 • Fax: 519-826-6947 • email: sales@IVFOonline.com