

# Voted #1 in IVF...

The Coda® Canister is the only filter with a 97% efficiency rating.  
With more activated carbon and potassium permanganate.

Homogenous blend of pure activated carbon and potassium permanganate removes odor, volatile organic compounds and chemical air contaminants.

Non toxic, no offgassing composite canister.

**Bad Air in...  
Good Air Out...  
Everytime...**

**Major reductions in:**

- Toxins
- Solvents
- Volatile Organics
- Chemical Contaminants
- Odororous Gases

**Removes more:**

- Toluene
- Benzene
- Hexane
- Formaldehyde
- Acetaldehyde
- Ethyl benzene
- Ethanol
- Isoproponal
- Acetone
- and more...

**Pure Clean Air**

**Pure Clean Air**

**IN**

Density and composition of media increases contact time and absorption.

Dual airtight seal ensures against leakage for both single and double pass systems.

CE

FDA 510(k) Cleared

The Coda® Canister has more absorption capacity and a longer life efficiency over all similar carbon, potassium permanganate filled filters in the market today. The IVF laboratory requires the added protection of the Coda® Canister in your Coda Tower®, Coda® Lo-Boy and HVAC systems.

*The IVF Laboratory requires Better Air Quality...*

*By genX the "Leader in Air Purification" in IVF.*




**IVFonline, LLC**  
Leading Fertility Brands

# The Coda® Canister Main Stage Filter

The Coda® Canister Main Stage Filter begins with a >95% absorption efficiency rate which is far above the 57% efficiency rate of conventional loose carbon filled V-cell filters.

Coda® Canisters, installing 5 in our custom frame, will adsorb almost 1 lb of VOCs and CACs, 2 times that of a V-cell filter while still maintaining the same efficiency rate.



The Coda® Canister will continually remove more harmful contaminants from the air in your IVF laboratory, without a reduction in air exchanges, protecting you and your specimens.

**The Coda® Canister easily loads into the Coda Tower®, Coda® Lo-Boy and most HVAC systems.**

Item	Frequency	Catalogue #
Coda® Canisters (5)	Six Months	CCST-005
Coda® Canister Frame	Reusable*	CCST-000
Coda® Canister Initial Kit Consists of: <ul style="list-style-type: none"> <li>• 5 Coda® Canisters</li> <li>• 1 Coda® Canister Frame</li> <li>• 1 Coda® HEPA Filter</li> <li>• 2 Coda® PreFilters</li> </ul>	Initial	For Coda Tower®: CCTI-081
		For Coda® Lo-Boy: CCLI-082
Coda® Canister Semi-Annual Kit Consists of: <ul style="list-style-type: none"> <li>• 5 Coda® Canisters</li> <li>• 2 Coda® PreFilters</li> </ul>	Semi-Annual	For Coda Tower®: CCTS-081
		For Coda® Lo-Boy: CCLS-082
Coda® Canister Annual Kit Consists of: <ul style="list-style-type: none"> <li>• 5 Coda® Canisters</li> <li>• 1 Coda® HEPA Filter</li> <li>• 2 Coda® PreFilters</li> </ul>	Annual	For Coda Tower®: CCTA-081
		For Coda® Lo-Boy: CCLA-082

\*The frame is reusable and you should order one at the initial purchase of the Coda® Canister System.

*Worldwide Service at Your Fingertips*

**Call or Click**



**US/Canada: 800-720-6375  
International: 519-826-5800**

**www.IVFOonline.com**

US/Canada: 800-720-6375 • International: 519-826-5800 • Fax: 519-826-6947 • email: sales@IVFOonline.com