



Maintaining Your Coda® Systems

Changing Your Filters Regularly Is Extremely Important



Coda Tower®
GTTU-040



V-Cell Quarterly Kits

GTFQ-041Coda Tower® Quarterly Kit
GXTQ-042Coda® Lo-Boy Quarterly Kit

Canister Semi-Annual Kits

CCTS-081Coda Tower® Semi-Annual Kit
CCLS-082Coda® Lo-Boy Semi-Annual Kit

**V-Cell Quarterly Kits – Please change every three months.
Canister Semi-Annual Kits – Please change every six months.**



Coda® Lo-Boy
GXTU-042



V-Cell Annual Kits

GTFA-041Coda Tower® Annual Kit
GXTA-042Coda® Lo-Boy Annual Kit

Canister Annual Kits

CCTA-081Coda Tower® Canister Annual Kit
CCLA-082Coda® Lo-Boy Canister Annual Kit

Annual Kits – Please change once per year.



Coda® Aero
GCAU-010



Coda® Aero Filters

GCAF-001Coda® Aero Filter (1 pack)
GCAF-002Coda® Aero Filter (2 pack)

(The Coda® Aero comes with an initial filter that must be changed after three months.)

Coda® Aero Filters – Please change every three months.



Coda® Unit
GCIU-010



Coda® Filter

GCIF-001Pkg of 1
GCIF-006Pkg of 6
GCIF-012Pkg of 12

Coda® SP
GCSP-010

Coda® Filters – Please change every thirty days.



Coda® Inline® Filter

GILF-001Pkg of 1
GILF-006Pkg of 6
GILF-012Pkg of 12

Coda® Xtra Inlines Filters®

GILF-001Pkg of 1
GILF-006Pkg of 6
GILF-012Pkg of 12

**Coda® Inline® Filters – Please change every two months or when cylinders are replaced.
Coda® Xtra Inline® Filters – Please change every three months or when cylinders are replaced.**

ISO9001:2000 and ISO13485:2003 Company

US/Canada: 800-720-6375
International: 519-826-5800
Fax: 519-826-6947
email: sales@IVFonline.com

IVFonline, LLC
innovation & quality

<genX>
international