

# Pioneer IVF

## Pioneer IVF gives you the 3 most reliable Aspiration Pumps for IVF...

Pioneer has provided the most reliable Pro-pump series available in IVF, a single vacuum, dual vacuum and the TriVac model. Pioneer Pro-pumps are used in a majority of IVF laboratories worldwide, available in both 115V and 230V.



### Single Vac Pioneer Pro-Pump

This standard aspiration pump is a powerful, durable Pioneer Pro-pump, operated with a **single** foot pedal. A fully adjustable vacuum control knob provides an even vacuum flow from 0 to 450 mmHg. Switch to 'high' with the push of a button on the front panel.

Pro-Pump Single (115V)	...GPPS-010/115
Pro-Pump Single (230V)	...GPPS-010/230
Accessory Kit	.....GPPK-075

### Dual Vac Pioneer Pro-Pump

Dual vacuum model lets you set a **'low'** vacuum level and to switch to a **'high'** vacuum simply by pressing specific foot pedal. A fully adjustable vacuum control knob provides an even vacuum flow from 0 to 450 mmHg.

Pro-Pump Dual (115V)	.....GPPD-010/115
Pro-Pump Dual (230V)	.....GPPD-010/230
Accessory Kit	.....GPPK-075



### TriVac Pro-Pump

Lets you set a **Level 1, Safe 'low' level** and a second **Level 2, Safe 'mid' range**, which is activated by the foot pedal. Lets you set **Level 1** at a very safe 40-80 mmHg vacuum and a **Level 2** level of 80-100 safe level with your foot pedal. Gives you a **Level 3 'high'** vacuum with the push of a button on the front panel.

Pioneer TriVac (115V)	.....PTPP-010/115
Pioneer TriVac (230V)	.....PTPP-010/230
Accessory Kit	.....PTAK-010

Patent Pending



Pioneer Pro-Pumps are the Most Reliable Aspiration Pump Used in IVF

- Quietest in the industry
- Light Weight
- Vacuum 0 to 300 mmHg
- Responsive
- 5 Year Warranty
- Made in USA
- Steel Housing
- 115V & 230V
- CE Register
- ISO 9001:2000
- ISO 13485:2003
- FDA 510(k) Cleared

# Pioneer IVF



[www.IVFOonline.com](http://www.IVFOonline.com)

US/Canada: 800-720-6375 • International: 519-826-5800 • Fax: 519-826-6947 • email: sales@IVFOonline.com